



# Flo-Form® Fluidised Positioners



## What is Flo-Form?

With its fluidised technology, Flo-Form can be moulded to any shape and will hold its form until it is remoulded. With its unique characteristic, Flo-Form is an excellent medium for patient positioning and pressure redistribution. Fluidised positioning devices (FPD's) are comprised of lubricated microspheres, giving them the ability to 3 - dimensionally contour to the patients specific shape



- Supports the neck and shoulders to offload the occiput
- Supports the Cervical Spine
- Promotes neutral Spine alignment
- Facilitates tube/line management
- Protects the ears, nose and face against shear and friction

## For Head Positioning



## For Body Positioning

- Contours and supports the body to offload the Scapula and Sacrum
- Maintains the Therapeutic Positioning
- Promotes neutral body alignment
- Facilitates tube/line management
- Moulds to patients body through 3D micro-contouring

